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**Purpose of the device:**

Static high resilience foam mattress intended to be used in conjunction with an adult medical bed for specific medical purposes due to a loss of autonomy and/or deterioration in health status that reduces mobility and exposes to the risk of pressure ulcers.

**Indications:**

This pressure ulcer prevention mattress is intended for patients with a medium to high risk of pressure ulcers (according to the Braden scale or another validated scale, and according to the doctor's advice).

Mattress used as a bed support to assist in the treatment of pressure ulcers or in the post-surgery treatment of pressure ulcers:

- In patients with one or more stage 1 and/or 2 pressure ulcers outside the support zone, or with the possibility of exclusion of support, whether the patient gets up during the day,
- In patients with a stage 3 or 4 pressure ulcer outside the support zone, or with the possibility of exclusion of support with a localised discharge system,
- In patients with a stage 1 or 2 pressure ulcer in the support zone, the patient is lifted during the day with a positioning system and medical auxiliaries are called in 3 times a day to check the installation and turn the patient over.

**Target patient population:**

Patients awake during the day and bedridden for more than 15 hours per 24 hours.

**Contraindications:**

- Patients who cannot be alternately positioned,
- Patients at very high risk of pressure ulcers,
- Patients with pre-existing pressure ulcers,
- Patients with a maximum weight greater than that tolerated by the mattress.

**Indicate undesirable side effects:**

Any serious incident occurring in connection with the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established. Inform the competent authority if you consider or have reason to believe that the device presents a serious risk or is a falsified device.

**Composition of the medical device:**

References	04.515D, 04.515DP	04.515DD, 04.515DDP	04.515T	04.515TD
Description	UNIMAT – BED 90		UNIMAT – 4 TRUNCATED CORNERS – BED 90	
Cover	CARTEX <sup>1</sup>	PHARMATEX bielastic <sup>2</sup>	CARTEX <sup>1</sup>	PHARMATEX bielastic <sup>2</sup>
Dimensions (cm)	88 x 198 x 15 cm			
Patient weight (kg)	Minimum and maximum: 30 to 140 kg			
Foam specification	Top layer	High resilience polyurethane foam (HR) 40 kg/m <sup>3</sup> (1,6 kPa)		
	Bottom layer	High resilience polyurethane foam (HR) 37 kg/m <sup>3</sup> (3,5 kPa)		
	Insert	High resilience polyurethane foam (HR) 40 kg/m <sup>3</sup> (6 kPa)		
Fire standards	NF EN 597-1 & NF EN 597-2			

<sup>1</sup>: 65% polyurethane, 35% polyester (weight 140 g/cm<sup>2</sup>); <sup>2</sup>: 50% polyurethane, 50% polyester (weight 205 g/cm<sup>2</sup>).

**Clinical benefits, performance, mechanism of action:**

*Performance characteristics of the device:* The foam mattresses allow for pressure relief in supporting areas to improve blood circulation and thus reduce the risk of pressure ulcers in the weight bearing areas of the body, through immersion and envelopment.

*Expected clinical benefits:* Maintenance of tissue oxygenation in the anatomical areas which are in contact with the supporting surface, by reducing the pressure applied to the cutaneous and sub-cutaneous tissue.

*Information for health professionals:* Observe the condition of the patient's skin which is in contact with the support structure several times a day. Use discharge devices or positioning systems on patients with pressure ulcers.

**Pre-use requirements and instructions for use:**

*User training:* User training must be carried out by people trained and approved by the economic operators concerned, particularly in terms of safety and reporting non-conformities.

*Installation of the device:* The mattress is delivered with its zipped cover. The product is ready to be installed. The cover must be dry. The head/feet symmetry principle enables the mattress to be positioned in 2 different positions. The cover is marked on the top, enabling teams to regularly alternate the direction of the mattress to enhance patient comfort.

*Preventive maintenance:* Carry out regular visual checks on the condition of the foam: visible sagging of the material and

uneven, very slow return of the foam are ageing criteria that compromise the properties of the substrate. Check the condition of the cover, the foam and the surface appearance of the cover annually. Annually check the outer surface of the cover by exposing the inside to a light source to check for holes and/or tears.

**Warnings, precautions, and required measures:**

*Precautions for use*

Non-stabilised bone injuries and/or muscular injuries in contact with the support. First days of post-surgical pressure ulcer treatment (skin graft or flap). Patients cared for at home with no possibility of medical auxiliary intervention. Also check the condition of exposed skin at each treatment and change of position. COPD-type bronchopulmonary disorders with a severe reduction in respiratory muscle tone. Joint stiffening and tendon retractions. Fractures of the pelvis or spine.

*Warnings*

Re-evaluate the patient’s risk of entrapment in the non-moving parts of the medical bed associated with “therapeutic mattress and accessories and articulated positions of the bed frame” in accordance with the EN 60601-2-52 standard for adults. The use of segmented or full-length side rails must be prescribed.

Non-compatible external cardiac massage without a board between the chest and the upper surface of the mattress.

The use of dry hot water bottles heated in the microwave in direct contact with the surface of the mattress is prohibited.

*Required measures*

- A mattress must be installed on a good quality flat base.
- Ensure that the dimensions of the mattress are compatible with the bed. The device must be used with its original protection.
- If necessary, check that the height of the bed rails is compatible with the thickness of the mattress (a minimum distance of 22 cm between the top of the bed rail and the top of the mattress when the patient is present is recommended).

Please note that support alone is not enough to prevent pressure ulcers/injuries. Other preventive measures are also essential, and are as follows:

- Change position at least every 2 to 3 hours,
- Maintain skin hygiene and avoid maceration. Do not tuck the bed linens too tightly; use elastic sheets to preserve mattress effectiveness,
- For the support to be effective, it is important to minimize the thickness between the body and the support as much as possible, except for the bedsheet in the case of a bed support, the undergarment, and a potential complete change. Prefer loosely fitted cotton undergarments and, if possible, those without seams in pressure areas. Do not insert folded towels or sheets, additional pillows, etc.,
- Alternate with static or dynamic standing, if possible,
- In case of incontinence, regularly change the pads,
- Observe or have observed daily the state of the skin,
- Ensure that the diet is sufficient and appropriate,
- Drink regularly and in sufficient quantity.
- If any of these measures cannot be followed, it is essential to inform your doctor or nurse as soon as possible.

*Circumstances in which the user should consult a healthcare professional*

Report to your doctor or nurse as soon as possible any abnormal event such as fever, pain, redness or whitening of the support points (head, shoulder, back, hip, shoulder blade, pelvis, heel, etc.).

*Information regarding suitable procedures for re-use*

The medical device may be reused as many times as possible, as long as its performance is still guaranteed and optimal.

*Accessories*

Only the mattress protection accessory supplied by the manufacturer ensures the safety and performance of the complete product.

**Cleaning and disinfection:**

Cover					
	Moderate wash up to 90°C.	Maximum allowable chlorine concentration of 5000 ppm.	Trumble dry on low heat.	Do not iron.	Do not dry clean.
Mattress	Cleaning with neutral detergent/disinfection with suitable sprays or aerosols.				

**Storage, handling, disposal:**

*Conditions of use and storage*

	Use	Storage
<b>Temperature range</b>	+15°C à +45°C / +59°F à + 113°F	-25°C à +70°C / +13°F à +158°F
<b>Humidity range</b>	30% - 70%	30% - 95%
<b>Atmospheric pressure</b>	50 kPA – 106 kPA	

Mattress should preferably be stored flat, away from direct light and excessive humidity.

*Product disposal:* Do not throw the product out into non-dedicated locations in nature. Use the recycling facilities available in your country.

*Lifetime:* The estimated lifetime of the product is 6 years.

Siotec AG, Sägestrasse 75, CH-3098 Köniz  
CHRN-AR-20001368