

IN GEL AND FOAM **ALOVA™**



**GELSCAR MINI HEEL HOLDER**



- ▶ Relieves pressure applied on heel
- ▶ Reduces the risk of shears
- ▶ Minimizes the risks of heating of the skin by friction

**SINGLE COMPARTMENT MINI HEEL HOLDER**



- ▶ Total heel support
- ▶ Ergonomic (easy to use)
- ▶ Compact (takes very little space)
- ▶ "Single patient" use


**FOOT LIFT SPLINT**



- ▶ Help in preventing heel pressure ulcers or in treatment of existing heel pressure ulcer
- ▶ Full heel relief
- ▶ Orthopedic alignment of the lower limb
- ▶ Lateral support and bending of the foot to prevent equinus
- ▶ Slight tilt of lower limb to encourage venous return



**STABILIZING HEEL HOLDER**



- ▶ Help in preventing heel pressure ulcers or in treatment of existing heel pressure ulcer
- ▶ Full heel relief
- ▶ Slight tilt of lower limb to encourage venous return



**WIDE HEEL HOLDER**



- ▶ Help in preventing heel pressure ulcers or in treatment of existing heel pressure ulcer



**HEEL HOLDER FOR BED BOTTOM**



- ▶ Relieves the heel area by respecting the physiological alignment of the lower limb.

**LATERAL DECUBITUS CUSHION**



- ▶ Allows positioning the patient in semi-lateral position at 30° to preserve areas at risk (sacrum and trochanter) by reporting the pressure on areas at lower risk of pressure ulcers without bony prominences and well vascularized. In addition, heel relief can be used for patient at high risk of pressure ulcers.

TECHNICAL AID DEVICES IN GEL AND/OR FOAM ALLOW PREVENTING FROM HEEL PRESSURE ULCERS BY RELIEVING PRESSURE POINTS AND VICIOUS POSTURES OF LOWER LIMBS IN LATERAL DECUBITUS POSITION

**KEY BENEFITS**

- ▶ Fit with different morphologies
- ▶ Allow relieving pressure points in area(s) at risk

DESCRIPTION	REFERENCE	LENGTH LEG SEGMENT	MATERIALS	PROTECTORS	USE
GELSCAR MINI HEEL HOLDER	VGLS13U	Heel location	FLUID GEL	—	One size
MINI HEEL HOLDER	VTM/54	19 cm/7.4 inch	ZYPREX	Without protection	
FOOT LIFT SPLINT	VTAC410/CIC	41cm/16.1 inch	ALOVA	PROMUST CIC	Reusable
	VTAC450/CIC	45 cm/17.7 inch			
STABILIZING HEEL HOLDER	VDT41A/CIC	41cm/16.1 inch			
	VDT45A/CIC	45 cm/17.7 inch			
WIDE HEEL HOLDER	VTLA/HIPH	45 cm/17.7 inch		Promust PU HD	
HEEL HOLDER FOR BED BOTTOM	VTLAXL/HIPH	45 cm/17.7 inch			
LATERAL DECUBITUS CUSHION	VCDA/50			Promust PU HD	

2 years warranty

Maintenance: Promust PU HD covers



Photos not contractually binding

FP000046/01\_2018



WINNOCARE Group - 200 rue Charles Tellier - Actiparc de Grézan - 30034 Nîmes Cedex 1 France  
 Ph : +33 (0)4 66 02 15 15 - Fax : +33 (0)4 66 02 15 00 - Email : contact.askle@winnocare.fr - www.winnocare.fr  
 In the interest of constantly improving its products, the WINNOCARE Group reserves the rights to modify the technical specifications of any product without prior notice  
 Class I medical device in accordance with European Directive 93/42/EEC

